- |           |   |  |           |  |  |
|-----------|---|--|-----------|--|--|
| <b>EN</b> | 1. Arm Cuff   | 5. Memory Key  | <b>ES</b> | 1. Brazalete   | 5. Botón de memoria  |
|           | 2. LCD Display  | 6. ON/OFF/START Key  |           | 2. Pantalla LCD  | 6. Botón de ENCENDIDO/ APAGADO/INICIO  |
|           | 3. Air Tube and Connector   | 7. Battery Cover   |           | 3. Tubo de aire y conector   | 7. Tapa de pilas   |
|           | 4. Inflating Bulb   |  |           | 4. Pera de inflado   |  |
| <b>DE</b> | 1. Oberarmmanschette, Luftschlauch zum Anschluss, Gummiball zum Aufpumpen | 5. „Speicher“-Taste, Taste ON/OFF/ START, Batterieabdeckung, Touchscreen | <b>FR</b> | 1. Brassard pour bras, écran LCD, Tuyau d'air et connecteur, Poire gonflable | 5. Touche Mémoire, Touche ON/OFF/ DEMARRAGE, Couverture du logement des piles, |

## English

you that you should continue pumping the pressure to more than 200 mmHg.

3. Sit still and wait for the monitor to deflate and measure your blood pressure. It is important to remain still and quiet during measurement. Any significant movement may affect measurement results.

4. When the measurement is completed, systolic blood pressure, diastolic blood pressure, and pulse will be shown simultaneously and be saved automatically in the memory system. Up to 90 measurements can be saved in the memory system.

5. Press the ON/OFF key to rapidly deflate and turn off the power. If no key is pressed, the unit will shut off automatically in 1 minute.

Note: This monitor automatically switches off approximately 1 minute after last key operation. To interrupt the measurement, simply press the Memory key or ON/OFF key; the cuff will deflate immediately. During the measurement do not talk or move your arm or hand muscles.

**Recalling Values from Memory**  
To read memory values, press the Memory key. The memory values will show on the LCD display. The last measurement will be shown first. Every new press of the Memory key calls for one prior value stored. When the Memory key can store up to 90 readings. When the memory capacity exceeds 90 readings, the oldest data will be replaced with the new record.

**Clearing Values from Memory**  
Press and hold the Memory key for approximately 5 seconds, then the data in the memory zone can be erased automatically.

**Troubleshooting**  
If any abnormality arises during use, please check the following points.

Symptoms	Check Points	Correction
No display when the ON/OFF key is pressed	Have the batteries run down? Have the batteries' polarity been positioned incorrectly? Is the cuff placed correctly?	Replace them with four new batteries. Re-check the batteries in the correct positions. So that it is positioned correctly.
EE mark shown on display	Is the cuff pressure pumped enough? (50-60 mmHg above normal systolic value)? Did you talk or move during measurement? Did you vigorously shake the cuff during measurement?	Measure it again and pump up 50-60 mmHg more than normal systolic value. Keep it quiet. Measure it again. Keep arm steady during measurement.

Note: If the unit does not work, return it to your dealer. Under no circumstance should you disassemble and repair the unit by yourself.

**Cautionary Notes**

- The unit contains high-precision assemblies. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or strongly shaking the main unit, and protect it from dust.
- Clean the blood pressure monitor body and the cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash the cuff or use chemical cleaner on it. Never use alcohol or petrol (gasoline) as cleaner.
- Leaky batteries can damage the unit. Remove the batteries when the unit is not used for a long time.
- The unit should not be operated by children so to avoid hazardous situations.
- If the unit is stored near freezing, allow it to acclimate at room temperature before use.
- This unit is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**Display Explanations**  
**EE / Measurement Error:** Make sure the L-Plus is securely connected to the inflating bulb and measure again quietly. Wrap the cuff correctly and keep arm steady during measurement. If the error keeps occurring, return the device to your local distributor or service center.

**E1 / Air Circuit Abnormality:** Make sure the L-Plus is securely connected to the air socket on the side of the unit and measure again quietly. If the error still occurs, return the device to your local distributor or service center for help.

**E2 / Pressure Exceeding 300 mmHg:** Switch the unit off and measure again quietly. If the error keeps occurring, return the device to your local distributor or service center.

**1.** Press down on the battery cover in the direction of the arrow to open the battery compartment.  
**2.** Install or replace 4 "AAA" sized batteries in the battery compartment according to the indications inside the compartment.  
**3.** Replace the battery cover by clicking in the bottom hooks first, then push in the top end of the battery cover.  
**4.** Replace the batteries in pairs. Remove batteries when unit is not in use for extended periods of time.

**1.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**2.** As a common issue for all blood pressure monitors using the oscillometric measurement function, the device may have difficulty in determining the proper blood pressure for users diagnosed with common arrhythmia (atrial or ventricular premature beats or atrial fibrillation), diabetes, poor circulation of blood, kidney problems, or for users suffered from stroke, or for unconscious users.

**3.** This unit is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**4.** Replace the batteries in pairs. Remove batteries when unit is not in use for extended periods of time.

**5.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**6.** This unit is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**7.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**8.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**9.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**10.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**11.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**12.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**13.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**14.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**15.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**16.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**17.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**18.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**19.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**20.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**21.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**22.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**23.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**24.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**25.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**26.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**27.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**28.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**29.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**30.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**31.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**32.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**33.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**34.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**35.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**36.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**37.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**38.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**39.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**40.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**41.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**42.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**43.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**44.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

**45.** Do not use this device for diagnosis or treatment of any health problem. This device is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact the store or the doctor from whom you purchased this unit or please contact Rossmax International Ltd.

\*Specifications are subject to change without notice.

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