



# rossmax

just a heartbeat away

prevention • monitoring • therapy • professional

## Dual Head Stethoscope

Model : EB200



- Soft ear tips
- Aluminium chest piece
- Non-chill bell and diaphragm ring
- Diaphragm side, detects high frequency sound
- Bell side, detects low frequency sound
- Snap on diaphragm ring
- Latex-free tube



Model	Qty per carton	Carton volume
EB200	60 pcs	0.052 cbm/ctn

www.rossmax.com  
Epaper\_EB200\_1\_ver1541